



ICON SELF-LIGATING BRACKETS

STD TORQUE

MAXILLARY	TORQUE	TIP	R/L	COLOR ID	HOOK	.022
Central (U1)	+12°	+5°	Right	●	-	3034-22-11
			Left	●	-	3034-22-21
Lateral (U2)	+8°	+9°	Right	●	-	3034-22-12
			Left	●	-	3034-22-22
Cuspid (U3)	+7°	+8°	Right	●	hook	3034-22-13H
			Left	●	hook	3034-22-23H
1st Bicuspid (U4)	-8°	+0°	Right	●	-	3034-22-14
			Left	●	-	3034-22-24
	-8°	+0°	Right	●	hook	3034-22-14H
			Left	●	hook	3034-22-24 H
2nd Bicuspid (U5)	-8°	+0°	Right	●	-	3034-22-15
			Left	●	-	3034-22-25
	-8°	+0°	Right	●	hook	3034-22-15H
			Left	●	hook	3034-22-25H

STD TORQUE

MANDIBULAR	TORQUE	TIP	R/L	COLOR ID	HOOK	.022
Central (L1)	-2°	+0°	Right	●	-	3034-22-41
			Left	●	-	3034-22-31
Lateral (L2)	-2°	+0°	Right	●	-	3034-22-42
			Left	●	-	3034-22-32
Cuspid (L3)	+7°	+3°	Right	●	hook	3034-22-43H
			Left	●	hook	3034-22-33H
1st Bicuspid (L4)	-12°	+0°	Right	○	-	3034-22-44
			Left	○	-	3034-22-34
	-12°	+0°	Right	○	hook	3034-22-44H
			Left	○	hook	3034-22-34 H
2nd Bicuspid (L5)	-17°	+0°	Right	●	-	3034-22-45
			Left	●	-	3034-22-35
	-17°	+0°	Right	●	hook	3034-22-45H
			Left	●	hook	3034-22-35H

FULL CASE

.022

5-5 hooks with hooks on 3

3034-22-553

5-5 hooks with hooks on 3,4,5

3034-22-55345